




**At the heart of life.
At the heart of precision.**






MyVal-1 Study 6-month Outcome

Low^{*} Device Related Mortality
Incidence of Stroke
New Permanent Pacemaker

97% Device Success^{*}

Myval THV: Designed for Precision in Outcomes





Myval THV has been indigenously developed by Meril Life Sciences Pvt. Ltd.

Myval THV: Unique Crimping Outcome



Upon Crimping:


- V-shaped hinges on hexagonal frame fold, generating the dense bands on fluoroscopy
- Vertical connectors 'I' give rise to light bands
- Alternating V-folds & vertical connectors give Myval THV a unique '**Zebra Crossing**' appearance on fluoroscopy for ease of positioning



Myval THV is recommended to be crimped over Navigator THV Balloon Delivery System prior to insertion within patient's vasculature.

Myval THV: Precise Placement Technique

Schematic of Myval THV - Ideal Landing Zone






Elimination of THV frame parallax promptly ensures visualisation of characteristic dark-light bands

Myval THV: Ground Zero Deployment

- Shallow deployment of Myval THV with least engagement within LVOT is possible.
- Optimal orthotopic anchorage of Myval THV with marginal LVOT foot-print without risk of THV migration.
- Minimal infra-annular depth $\leq 3.50\text{mm}$ avoids conduction system interference (thus minimizing the need of new permanent pacemaker dependency).

Largest circumscribable diameter in Open Cell
(for all Myval THV Diameters 20mm to 32mm)







Myval THV infra-annular depth $\leq 3.50\text{mm}$ based on THV \varnothing

Myval THV: Detailed Sizing Guide

3D Annular area mm ²		270	280	290	300	310	314	320	330	340	350	363	
3D area derived diameter mm		18.5	18.9	19.2	19.5	19.9	20.0	20.2	20.5	20.8	21.1	21.5	
% Annular area over/under	20 mm	16.4%	12.2%	8.3%	4.7%	1.3%	0.1%	-1.8%	-4.8%	-7.6%	-10.2%	-13.5%	
	21.5 mm	34.5%	29.7%	25.2%	21.0%	17.1%	16%	13%	10%	7%	4%	0.0%	
	23 mm	53.9%	48.4%	43.3%	38.5%	34%	32.3%	29.8%	25.9%	22.2%	18.7%	14.5%	
3D Annular area mm ²		370	380	390	400	410	415	420	430	440	450	460	471
3D area derived diameter mm		21.7	22.0	22.3	22.6	22.8	23.0	23.1	23.4	23.7	23.9	24.2	24.5
% Annular area over/under	23 mm	12.3%	9.3%	6.5%	3.9%	1.3%	0.1%	-1.1%	-3.4%	-5.6%	-7.7%	-9.7%	-11.8%
	24.5 mm	27.4%	24.1%	20.9%	17.9%	15.0%	13.6%	12.2%	9.6%	7.1%	4.8%	2.5%	0.1%
	26 mm	43.5%	39.7%	36.1%	32.7%	29.5%	27.9%	26.4%	23.5%	20.7%	18.0%	15.4%	12.7%
3D Annular area mm ²		480	490	500	510	520	531	540	550	560	570	580	594
3D area derived diameter mm		24.7	25.0	25.2	25.5	25.7	26.0	26.2	26.5	26.7	26.9	27.2	27.5
% Annular area over/under	26 mm	10.6%	8.4%	6.2%	4.1%	2.1%	0.0%	-1.7%	-3.5%	-5.2%	-6.9%	-8.5%	-10.6%
	27.5 mm	23.7%	21.2%	18.8%	16.5%	14.2%	11.9%	10.0%	8.0%	6.1%	4.2%	2.4%	0.0%
	29 mm	37.6%	34.8%	32.1%	29.5%	27.0%	24.4%	22.3%	20.1%	17.9%	15.9%	13.9%	11.2%
3D Annular area mm ²		600	610	620	630	640	650	661	670	680	690	700	710
3D area derived diameter mm		27.6	27.9	28.1	28.3	28.5	28.8	29.0	29.2	29.4	29.6	29.9	30.1
% Annular area over/under	29 mm	10.1%	8.3%	6.5%	4.8%	3.2%	1.6%	-0.1%	-1.4%	-2.9%	-4.3%	-5.6%	-7.0%
	30.5 mm	37.6%	37.6%	37.6%	37.6%	14.2%	12.4%	10.5%	9.0%	7.4%	5.9%	4.4%	2.9%
	32 mm	34.0%	31.8%	29.7%	27.7%	27.7%	23.7%	21.7%	20.0%	18.3%	16.6%	14.9%	13.3%
3D Annular area mm ²		720	731	740	750	760	770	780	790	804			
3D area derived diameter mm		30.3	30.5	30.7	30.9	31.1	31.3	31.5	31.7	32.0			
% Annular area over/under	32 mm	11.7%	10.0%	8.7%	7.2%	5.8%	4.4%	3.1%	1.8%	0.0%			

Area derived diameter based on MSCT cross-sectional measurements.

Myval THV: Size Matrix


Myval THV Size Matrix & Technical Specifications	Area 314 mm ² 17.35 mm  20 mm	Area 363 mm ² 18.35 mm  21.5 mm	Area 415 mm ² 17.85 mm  23 mm	Area 471 mm ² 18.75 mm  24.5 mm
Perimeter	62.83 mm	67.54 mm	72.26 mm	76.97 mm
Native annulus area	270 - 330 mm ²	314-380 mm ²	360 - 440 mm ²	410-500 mm ²
Area-derived diameter	18.5 - 20.5 mm	20-22 mm	21.4 - 23.7 mm	22.8-25.2 mm
Native annulus size by TEE	16 - 19 mm	17.5-20.5 mm	18 - 22 mm	19.5-23.5 mm

All Myval THV diameters (20 mm to 32 mm) are compatible with 14Fr Python - Introducer Sheath

← Myval THV XL Sizes →

Area 531 mm ² 18.85 mm  26 mm	Area 594 mm ² 19.25 mm  27.5 mm	Area 661 mm ² 20.35 mm  29 mm	Area 731 mm ² 20.9 mm  30.5 mm	Area 804 mm ² 21.14 mm  32 mm
81.68 mm	86.39 mm	91.11 mm	95.82 mm	100.53 mm
460 - 560 mm ²	510-630 mm ²	570 - 700 mm ²	630-770 mm ²	700-840 mm ²
24.2 - 26.7 mm	25.5-28.3 mm	26.9 - 29.9 mm	28.3-31.3 mm	29.9-32.7 mm
21- 25 mm	22.5-26.5 mm	24 - 28 mm	25.5-29.5 mm	27-31 mm

Myval THV: Post Deployment Dimension Chart



Myval THV Diameters (Ø)	20 mm	21.5 mm	23 mm
Total frame height	17.35 mm	18.35 mm	17.85 mm
Open cell height (53%)	9.20 mm	9.73 mm	9.46 mm
Closed cell height (47%)	8.15 mm	8.62 mm	8.39 mm
Infra-annular depth	3.05 mm	3.20 mm	2.85 mm
Supra-annular height of closed cells	5.10 mm	5.42 mm	5.54 mm
Recommendation for coronary protection	10 mm	10 mm	10 mm

- A balloon occlusion test may be considered to assess the propensity for coronary occlusion. Balloon diameter approximated to shortest axis of CT derived annular diameter to be considered.

← Myval THV XL Sizes →

24.5 mm	26 mm	27.5 mm	29 mm	30.5 mm	32 mm
18.75 mm	18.85 mm	19.25 mm	20.35 mm	20.90 mm	21.14 mm
9.94 mm	9.99 mm	10.20 mm	10.79 mm	11.08 mm	11.21 mm
8.81 mm	8.86 mm	9.05 mm	9.56 mm	9.82 mm	9.94 mm
2.95 mm	3.05 mm	3.15 mm	3.35 mm	3.45 mm	3.55 mm
5.86 mm	5.81 mm	5.90 mm	6.21 mm	6.37 mm	6.39 mm
10 mm	10 mm	10 mm	10 mm	10 mm	10 mm

- Consider protection of coronary arteries with a DES especially if height of coronary ostium is < 10 mm from the annular plane and in conjunction with sinus of valsalva dimensions i.e. height & diameters.

Navigator THV Delivery System


Delivering TAVI Made Easy

- Myval THV is recommended to be crimped over Navigator THV Delivery System prior to insertion within patient's vasculature.
- The crimped valve with delivery system is then loaded through 14Fr Python – Introducer Sheath.



- Navigator delivery system has a set of proximal and distal stoppers which ensure that valve crimping is precise and snug.
- Visual confirmation of crimped valve can be ensured before entering the sheath to avoid any crimping errors/defects.
- The stoppers prevent inadvertent migration of the valve & ensure there is no risk of valve dislodgement (embolization) during entry through the sheath or while negotiating the loaded delivery system across the aorta.
- Myval THV direct crimping on the balloon makes TAVI delivery simple, intuitive and eliminates unwarranted procedural steps.

Navigator THV Delivery System



Navigator – THV Delivery System has been indigenously developed by Meril Life Sciences Pvt. Ltd.

Navigator THV Delivery System


Characteristic Balloon Expansion

Navigator balloon with dual expansion ports at each end ensures rapid, simultaneous, controlled expansion (dog-boning) of distal and proximal ends


This typical dog bone pattern of inflation steadies the valve during expansion phase, ensuring its precise annular position and deployment without any risk of valve migration

Rapid balloon inflation, using an inflation device is possible with controlled palm thrust


Rapid balloon deflation within 3-5 sec ensures procedural safety and compliance



Navigator THV Delivery System: Proximal Assembly



Hi-flexion feature ensures tracking the THV delivery system via inner aortic arch curve thereby avoiding contralateral wall scraping.



Caution: Always remember to fully un-flex the Navigator system while withdrawing


14Fr Python – Introducer Sheath

Compatible with all Myval THV diameters (20 mm to 32 mm)

Sheath expands momentarily like a python swallowing its prey
Conveniently allows passage of crimped Myval THV System

14Fr Entry Profile, Allows Atraumatic Percutaneous Access

High convenience for full retrievability of an un-deployed Myval THV System





Two separate, calibrated loading tubes ensure temporary opening of hemostatic valves in proximal port allowing smooth passage of crimped Myval THV System

Common Femoral Artery* Ø (mm)	Myval THV Ø (mm)
≥ 5.50 mm	20 mm, 21.5 mm, 23 mm, 24.5 mm
≥ 6.00 mm	26 mm, 27.5 mm, 29 mm
≥ 6.50 mm	30.5 mm, 32 mm
*CFA Ø must be MSCT derived. Excluding circumferential Ca ²⁺	

Python - Introducer Sheath has been indigenously developed by Meril Life Sciences Pvt. Ltd.

Mammoth - OTW Balloon Catheter

It is not mandatory to pre-dilate the native annulus prior to Myval THV implantation. Operators may consider pre-dilatation based on anatomical considerations.



Mammoth Balloon Ø	14 mm	16 mm	18 mm	20 mm	23 mm	25 mm
Balloon length	↔ 40 mm ↔					
Balloon RBP	↔ 6 atm ↔					
Volume of 75:25 saline:contrast to achieve stated diameter	7 ml	10 ml	13 ml	16 ml	23 ml	25 ml

Mammoth – OTW Balloon Catheter has been indigenously developed by Meril Life Sciences Pvt. Ltd.



~4000 Patients

Myval THV: Global Clinical Program

MyVal-1¹ (n=100)

Meril Initiated
Single arm
Prospective FIH study
29 Indian sites
STS $\geq 4 - \leq 15$
Primary 30-d endpoint achieved
**6m f/up presented
@ EuroPCR 2020**

Myval-EU² (n=250)

Compare-TAVI³ (n=1,062)

LANDMARK RCT⁴ (n=768)

Myval-Global⁵ (n= 2,000)

Myval-China⁶ (n=125)

²Meril Initiated
Single arm

Retrospective EU Registry

20 EU sites

Real world STS

Data collection in progress

**30-day f/up est.
@ EuroPCR 2021**

³Investigator Initiated

RCT – 1:1
Myval : Sapien 3 THV

Prospective study

10 Nordic + EU sites

Real world STS

Enrollment has started

**1-year f/up est.
@ EuroPCR 2023**

⁴Meril Initiated

RCT – 2:1:1
Myval : Sapien : Evolut THV series

Prospective study

50+ EU/ANZ sites

Real world STS

EC/IRB approvals in progress

**First patient first visit est.
@ QIV 2020**

⁵Meril Initiated
Single arm

Prospective study

100 Global sites

Real world STS

Under planning

**First patient first visit est.
@ QIII 2021**

⁶Meril Initiated
Single arm

Prospective study

10 China sites

Real world STS

Pre-study activities initiated

**First patient first visit
est. @ QI 2022**


MyVal-1: Study Design

A prospective, multicentre, single-arm, open-label study of Myval THV in the treatment of severe symptomatic native aortic valve stenosis.

Total number of patients: 100

Device Sizes – 20, 21.5, 23, 24.5, 26, 27.5 and 29 mm

CLINICAL FOLLOW-UP



Dr. Samin Sharma - Chairman
New York, USA



Dr. Ashok Seth - Principal Investigator
New Delhi, India



Dr. Praveen Chandra - Co-ordinating PI
New Delhi, India



Dr. Ravinder Singh Rao - Co-PI
Jaipur, India




Dr. P. K. Goel - Scientific Advisor
Lucknow, India

Study Investigators: Sharma Samin¹, Chandra Praveen², Ashok Seth³, Rao Ravinder Singh⁴, Goel P. K.⁵, Bharadwaj Prashant, Sethi Rishi, Sengottuvelu G., Mahajan Ajaykumar, Jose John, Abhaichand Rajpal, Ajit Kumar V K, Manjunath C N, Mehtrotra Sanjay, Rao Suryaprakash, Chaurasia Amit Kumar, Bahl V K, Kaul Upendra, Jain R K, Gopalumurugan AB, Rath P C, Trehan Vijay, Vivek Kumar, Roy Sanjeeb, Mantri R R, Sharma S M, Kler T S, Nair R C, Mehta Ashwin 1. Chairman; 2. Principal Investigator; 3. Co-ordinating PI; 4. Co-PI; 5. Scientific Advisor.
*6-month outcome data presented by Dr. Ravinder Singh Rao at EuroPCR 2020.
MyVal-1: Study (CTRI/2016/11/007512).

MyVal-1: Baseline Characteristics

Patient History	
Average Age (years)	73 ± 7.49
Mean STS	5.12%
History of Coronary Artery Bypass Graft surgery	17%
History of Previous PCI	13%
History of previous Aortic Valvuloplasty	1%
Cerebral vascular disease	3%
Peripheral vascular disease	1%



- Intermediate sizes were introduced after 90% of enrollment completion.

MyVal-1: Clinical outcomes up to 6-month follow-up

Excellent clinical safety & efficacy

Events	Post-procedure	1-Month Follow-Up	6-Month Follow-Up
Survival	98%	97%	91%
All-cause mortality	2%	3%	9%
Stroke	1%	2%	2%
Acute renal failure	2%	2%	2%
Life-threatening or disabling bleeding	1%	1%	1%
Endocarditis	0%	0%	1%
Myocardial infarction	0%	0%	0%
Major vascular complications	1%	1%	1%
Minor vascular complications	2%	2%	2%
Repeat hospitalization	NA	8%	10%
New permanent pacemaker	2%*	2%	2%

*One patient had RBBB pre-procedure


MyVal-1: Echocardiographic Findings at 6-month Follow-up

Echocardiographic findings

Parameters	Baseline	Post Procedure	30-day FU	6-month FU
Effective orifice area, (cm ²)	0.6 ± 0.2	1.7 ± 0.3	1.7 ± 0.5	1.8±0.5
Mean aortic-valve gradient, (mmHg)	47.4 ± 8.8	8.0 ± 2.7	8.8 ± 2.5	10.5±2.6
Peak aortic-valve gradient, (mmHg)	71.7 ± 13.0	14.4 ± 2.4	15.7 ± 2.8	17.9±2.9
Trans-aortic velocity, (m/s)	4.5 ± 0.4	1.9 ± 0.4	1.8 ± 0.4	1.8±0.3
Mean LVEF, (%)	45.5 ± 11.5	47.8 ± 11.1	48.6 ± 8.9	48.8±8.0
Moderate or severe mitral regurgitation, (n)	2	0	0	0
Aortic regurgitation, (n)	-	0	0	0


Values are mean ±SD (n) or % (n); FU: follow-up, LVEF: left ventricular ejection fraction.

Sustained Low Mean Gradients Post-Procedure and
~1.8cm² Large EOA at 6-month Follow-up (p<0.0001)




MyVal-1: Marked improvement in Quality of Life (QoL) parameters

NYHA Classification




MyVal-1: Marked improvement in Quality of Life (QoL) parameters

Six-minute walk test




Kansas City Cardiomyopathy Questionnaire Score



MyVal-1: Study Conclusion

- In 100 intermediate and high-risk patients of MyVal-1 FiH study, Myval THV system demonstrated excellent clinical and hemodynamic outcomes at 6-month:
 - 91% survival & low incidence of stroke (2%)
 - Low 2% rate of new permanent pacemaker implantation post-procedure
 - High procedural success (97%) due to precise orthotopic valve positioning
 - Significant improvement
- In real world global experience of ≈ 1000 cases; Myval THV has been consistently demonstrating high procedural success and clinical performance
 - Unique hybrid honey-comb geometry for precise positioning and orthotopic deployment.
 - Preserve THV geometry & respect patient's anatomy; Intermediate \varnothing 21.5, 24.5, 27.5 mm & XL \varnothing 30.5, 32 mm
 - Direct THV crimping on Navigator balloon makes TAVI delivery simple, intuitive and eliminates unwarranted procedural steps.
 - Compatibility of novel 14Fr Python Introducer sheath for all Myval THV \varnothing s; with high convenience of full retrievability of an un-deployed Myval THV system

LANDMARK RCT – 50+ Sites EU+ANZ



Primary Endpoint – 30 Days


All cause mortality | All stroke | Life threatening bleeding | Vascular complications
Acute Kidney Injury | Paravalvular leak (PVL) | New permanent pacemakers

ECG/Echo Follow-up
Video Densitometry
Clinical Follow-up

Baseline | Post Procedure | 30 D | 1 Y | 3 Y | 5 Y
Post Procedure
Upto 10 years

LANDMARK RCT EC/IRB work initiated. First patient enrollment expected Q4, 2020.

Myval THV System and Components




Myval - THV



Navigator - THV Delivery System



Mammoth - OTW Balloon



Python - Introducer Sheath

Myval THV System and Components - Ordering Information

Myval - THV Ordering Information

Diameters	20.0 mm	21.5 mm	23.0 mm	24.5 mm	26.0 mm	27.5 mm	29.0 mm	30.5 mm	32.0 mm
Product code	MVL200	MVL215	MVL230	MVL245	MVL260	MVL275	MVL290	MVL305	MVL320

Navigator - THV Delivery System Ordering Information

Diameters	20.0 x 30 mm	21.5 x 30 mm	23.0 x 30 mm	24.5 x 30 mm	26.0 x 30 mm	27.5 x 35 mm	29.0 x 35 mm	30.5 x 35 mm	32.0 x 35 mm
Product code	NVT20030	NVT21530	NVT23030	NVT24530	NVT26030	NVT27535	NVT29035	NVT30535	NVT32035

Mammoth - OTW Balloon Ordering Information

Diameters	16.0 x 40 mm	18.0 x 40 mm	20.0 x 40 mm	23.0 x 40 mm	25.0 x 40 mm
Product code	MTV1640	MTV1840	MTV2040	MTV2340	MTV2540

Python - Introducer Sheath Ordering Information

Product code	PHT14
--------------	-------

Val-de-Crimp - Heart Valve Crimping Tool Ordering Information

Product code	VLDC
--------------	------

Myval™, Navigator™, Python™, Mammoth™ & Val-de-Crimp™ are registered trademarks of Meril Life Sciences Pvt. Ltd.

For the use of a Registered Medical Practitioner of a Hospital or Laboratory only.

Illustrations are artist's representation and should not be considered as engineering drawings or photographs.

Meril's Global Presence

Meril Tibbi Cihazlar

Meril Tibbi Cihazlar İmalat ve Ticaret A.Ş.
İçerenköy Mah. Çetinkaya Sok. Prestij Plaza No:28
Kat:4 Ataşehir, 34752
İstanbul / Turkey
T: +90 216 641 44 24
F: +90 216 641 44 25

Meril Medical LLC.

Nauchnyi Proezd 19,
Moscow, Russia – 117 246.
Office - +7 495 772 7643

Meril GmbH.

Bornheimer Strasse 135-137,
D-53119 Bonn.
Germany.
T: +49 228 7100 4000
F: +49 228 7100 4001

Meril, inc.

2436 Emrick Boulevard,
Bethlehem, PA - 18020
T: + 610 500 2080
F: + 610 317 1672

Meril South America

Doc Med LTDA Al. dos Tupiniquins,
1079 - Cep: 04077-003 - Moema.
Sao Paulo, Brazil.
T: +55 11 3624 5935
F: +55 11 3624 5936

EU Representative.

Obelis S.A.Bd, General Wahis 53,
1030, Brussels, Belgium.
T : +32 2 732 5954
F : +32 2 732 6003
E mail@obelis.net

Manufacturer

Meril Life Sciences Pvt. Ltd.

Survey No. 135/139, Bilakhia House,
Muktanand Marg, Chala,
Vapi 396 191.
Gujarat, India.
T: +91 260 305 2100
F: +91 260 305 2125

Meril China Co. Ltd.

Room 2301B, 23F, Lixin Plaza,
No 90, South Hubin Road,
Xiamen China
T: 0086-592-5368505
F: 0086-592-5368519

Meril Life Sciences Pvt. Ltd.

301, A - Wing, Business Square,
Chakala, Andheri Kurla Road,
Andheri East, Mumbai 400 093
T: + 91 22 39350700
F: + 91 22 39350777
E: askinfo@merillife.com

Meril SA Pty Ltd

102, 104, S101 and S102,
Boulevard West Office Park,
142 Western Service Road,
Erf 813 Woodmead Extension 17
Sandton, Johannesburg – 2191
South Africa
T: +27 11 465-2049
F: +27 86 471 7941

**At the heart of life.
At the heart of precision.**

